


THE SPEED OF LIGHT



CAUTION:

- **Never look directly into the beam without a diffusing screen between you and the laser. If viewed directly, the beam from even a low power laser can cause permanent damage to your vision.**
- **Do not touch the surfaces of the mirrors. The mirrors are aluminized-front-surface mirrors and do NOT have a protective overcoating. Touching, or any attempt to clean the mirror surfaces will result in damage to the mirrors.**

INTRODUCTION

This experiment is a modification of Foucault's method for the measurement of the speed of light using a rotating mirror. In 1862, Foucault measured $c = (298\,000 \pm 500)$ km/s, or an accuracy of $\pm 0.17\%$ using a baseline of only 20 m. In 1929 Michelson used a development of the rotating mirror arrangement over a baseline of 22 miles between Mt. Wilson and Mt. San Antonio in California to determine $c = (299\,796 \pm 4)$ km/s. Although the technique you use is similar, you can probably only expect to achieve accuracies of $\pm 5\%$. You will find as did Foucault, that there is a limit to the accuracy achievable in a technique dependent on the measurement of the displacement of a light beam. It is interesting, however, considering the magnitude of the velocity being measured, that one can quite easily measure c in one or two laboratory periods.

The other challenge of this experiment is the solving of the problem of how to set-up an apparatus systematically. If you work out your system of alignment carefully, you can perform the experiment easily. However, if you do not proceed systematically, you could spend a couple of lab periods without obtaining any results.

THEORY

A schematic diagram of the apparatus is given in Figures 1 and 2.




Figure 1




Figure 2

Light from a source S falls onto the rotating mirror M_r and is focussed by the lens L to form an image I on the surface of the mirror M_o . The purpose of the intermediate mirror M_i (not shown in Figure 2) is to create a longer path by allowing the light to double back in the room. The light from the image on M_o will be reflected and form a second image exactly coincident with the original object S since I and S are conjugate points of the lens L . As a result of the focussing properties of L , all the rays that pass through L and reach M_o will reflect and form an image coincident with the object S even though the rays will not generally retrace their paths exactly.

An angle β is defined as the angle between the plane of the mirror and some arbitrary direction as shown in Figure 2. It is important to note that the position of the final image on return to the source location S is independent of the angle β of the rotating mirror M_r , although the position of the image I on the face of mirror M_o will vary with angle β . Of course, this is true only if the rays do pass through lens L and strike the surface of M_o . If they don't, the light just won't make it back to the location of S .

What is described above is true when M_r is stationary at any reasonable angle β or rotating at some low value of angular velocity $\omega = \frac{d\beta}{dt}$. However, if ω is sufficiently large so that the mirror M_r has rotated through an angle $\delta\beta$ in the time the light takes to travel from M_r to M_o and back, the final image will be displaced sideways from the original source position S .




Figure 3

Referring to Figure 3, if the rotating mirror is initially at position **1** then a line from S to P is perpendicular to M_r and the source S produces an image in M_r at P . The light, after reflection from M_o heads to P on its return and forms an image at S . But if M_r is now at position **2**, the light is reflected at a larger angle to form an image at S' . If the mirror has rotated through an angle $\delta\beta$, using the laws of reflection, you should be able to show that the angle subtended at the mirror M_r by S and S' is $2 \times \delta\beta$. Thus, the displacement x of the final image, due to change in mirror angle from when the light leaves M_r to when it returns is just $x = \ell \times 2\delta\beta$, assuming $\delta\beta$ to be small, and where ℓ is the distance from M_r to S .

If the optical path length between M_r and M_o is d , then the light, travelling at a speed c , takes a time $\delta t = \frac{2d}{c}$ between reflections at M_r , so that if the mirror is rotating at an angular velocity ω ,

then $\delta\beta = \omega \times \delta t = 2\omega \frac{d}{c}$ and thus

$$x = \frac{4\ell d}{c} \omega \quad (1)$$

The determination of c thus reduces to the measurement of two fixed distances, a variable angular velocity, and the displacement of an imaged light beam.

THE EXPERIMENT

Although the method requires merely a small intense light source we have provided you with a laser light source. Lasers have the advantage of lots of intensity, all in the direction you want.

SETUP

The rotating mirror M_r and the fixed mirrors M_o and M_i are placed to give the maximum path length possible given the finite size of the space. The rotating mirror M_r can be rotated by hand by inserting the removable knob into the top of the motor assembly. As a safety feature, the removable knob is tethered to an elastic cord. Removing the knob before starting the motor prevents accidental damage to the motor shaft.

Preliminary Beam Alignment

With lens L removed from the light path, use the following steps as a guide to form the path

$$S \Rightarrow M_r \Rightarrow M_i \Rightarrow M_o \Rightarrow M_i \Rightarrow M_r \Rightarrow S$$

1. Adjust the laser S and its stand so that the beam strikes the centre of M_r and reflects back to S .
2. Using the cement blocks in the wall of the room as a guide, adjust M_r and the laser so that the beam rotates in a horizontal plane as M_r is rotated. When making corrections to the alignment, make only partial corrections with either M_r or S . Note that when M_r is adjusted, S must also be simultaneously adjusted so that the beam still strikes the centre of M_r and still reflects back to S . This procedure may require several iterations.
3. Check that the centre of M_o is at the same height as the centre of M_r . If it is not, make the necessary adjustments.
4. Adjust M_i so that the beam hits its centre and use the adjustment screws so that the reflected beam strikes the centre of M_o .
5. Adjust M_o so that the reflected beam hits the centre of M_i and, even though it is large, the centre of the beam reflecting from M_i again strikes M_r . You should see a large red spot which is a rectangular image of M_r back at the laser source S .

Introduction of the Lens L

6. Calculate the position for lens L to make S and M_o optically conjugate (i.e., L causes an object at S to form an image at M_o). Use the thin lens formula $\frac{1}{s} + \frac{1}{s'} = \frac{1}{f}$ where s and s' are the object and image distances respectively. You will need to know the total distance from S to M_o . The focal length f , of L , is given in millimetres on the lens holder. Although the correct answer requires the solution of a quadratic equation, you should be able to estimate the correct answer in your head.

- Place **L** in its calculated location and adjust its vertical and transverse position so that the beam now traverses the path

$$\mathbf{S} \Rightarrow \mathbf{M}_r \Rightarrow \mathbf{L} \Rightarrow \mathbf{M}_i \Rightarrow \mathbf{M}_o \Rightarrow \mathbf{M}_i \Rightarrow \mathbf{L} \Rightarrow \mathbf{M}_r \Rightarrow \mathbf{S}$$

Introduction of the "Scoop Off" Mirror \mathbf{M}_s

By this point in the set-up, you have the beam returning to the laser source **S**. If you were to now run the motor to rotate \mathbf{M}_r you would not be able to make measurements since this would require that your eye be inside the laser. Thus you must now insert the "scoop-off" mirror \mathbf{M}_s between the laser **S** and \mathbf{M}_r , angled in such a way that the returning light beam gets reflected on to your observing scale rather than returning to the laser **S**.

Since the beam is very intense, use the coloured filter provided to reduce the brightness when viewing the beam through the magnifying eyepiece. Note that the smallest divisions on the viewing scale are 0.5 mm. The coloured filter is to be removed later when \mathbf{M}_r is in motion.

- The "scoop-off" is achieved by first making a very slight adjustment of one foot of the laser stand so that the beam strikes \mathbf{M}_r near its top. Blocking \mathbf{M}_o will make this adjustment easier.
- Probably the returning beam will be below \mathbf{M}_r (there may be multiple reflected beams but the returning beam will be the brightest). Make a very slight adjustment of the tilt of \mathbf{M}_o so that the returning beam strikes \mathbf{M}_r near its bottom. Use a piece of paper near \mathbf{M}_r to view the separation between the outgoing and returning beams.
- Now \mathbf{M}_s may be placed near \mathbf{M}_r such that it intercepts only the beam travelling from \mathbf{M}_r headed towards **S**, but does not intercept any of the three other light beams arriving at or leaving \mathbf{M}_r . Adjust the angle of \mathbf{M}_s so that the outgoing beam hits the viewing scale.
- Position the viewing scale at such a distance that the light spot is a minimum size. Although you should not move \mathbf{M}_r now, it should be pointed out that at this position the spot will not move from side to side when \mathbf{M}_r is rotated through small angles. Later when \mathbf{M}_r is rotating, the beam hits the scale for only a small range of angles and this is why you see a circular spot and not a flattened ellipse.
- Turn on the autotransformer and set the motor rotating at a low speed. Remove the coloured filter from the beam path. Using the eyepiece, fine tune the position of the scale so that the spot on the scale is round and small in size.

Measurements

- Note that the definition of ℓ now becomes the distance from \mathbf{M}_r to \mathbf{M}_s to the viewing scale. When the viewing scale is placed at the position where the spot size is a minimum, theory says that the light path distance from \mathbf{M}_r to \mathbf{M}_s to the viewing scale is the same as the distance from **S** to \mathbf{M}_r . Measure ℓ .

13. Measure **d**.

14. The time of rotation of the mirror and hence ω can be measured using the photodiode. Connect the photodiode to **CH1** of the oscilloscope. Press **CH1 Menu** and set:

Coupling - DC

BW Limit - ON

VOLTS/DIV - Coarse

Probe – 1X

INVERT - OFF

With **M_r** stationary and with the detector in a position such that no part of the beam is intercepted by the motor housing, adjust the photodiode for maximum signal which should be greater than 5 volts.

15. To determine the frequency when **M_r** is rotating, press **MEASURE**, choose **Type** and set **CH1** to **Freq**. Adjust **VOLT/DIV**, **SEC/DIV** and **TRIGGER LEVEL** as the motor speed changes to keep at least two peaks on the screen and the frequency measurement valid.

Without touching the mirror surfaces, when you determine ω , be sure to check whether the rotating mirror is single or double sided as this will affect your calculations.

16. Measure **x** as a function of ω over as large a range of ω as is practical. Take some of your data with ω increasing and some with ω decreasing and keep track of the order in which data is recorded. At high frequencies the photodiode may have to be adjusted for maximum signal.

DATA INTERPRETATION

The form of equation (1) suggests a graphical way of interpreting the data. The **x** in equation (1) represents the displacement of the light spot from its position when the mirror is not rotating. The measuring scale will probably be arbitrarily positioned with its zero not corresponding to the position of the spot for a stationary **M_r**. Thus equation (1) is better written as:

$$\mathbf{x} = \frac{4\ell\mathbf{d}}{\mathbf{c}}\omega + \mathbf{x}_0 \quad (2)$$

Do not try to determine **x₀** with the mirror stationary as the brightness of the spot is so much greater than when the mirror is rotating that the data taken under the two intensity conditions are not comparable.

SUPPLEMENT

One may be able to increase the precision of the measurements by inserting a slit directly in front of the laser. The slit will now be the source. Choose the slit width and slit position to optimize clarity, sharpness and intensity of image as seen in the viewing eyepiece. Note that it will be more difficult to measure higher frequencies on the oscilloscope since the intensity of the light striking the photodiode will now be lower.

(jv – 1988, jp - 2006)